

Permobil

F5 Corpus[®]



permobil

Agile everywhere ability

The F5 Corpus® is ideal for an active outdoor lifestyle and for driving on varying surfaces or over longer distances on regular basis. With its tight turning circle it offers notable manoeuvrability.

Corpus® – follows the contours of the body

The Corpus® seating system features an ergonomic design that follows the contours of the body. Improved comfort, pelvic positioning, lateral stability and immersion are all key objectives of the new Corpus cushions.

Corpus cushions feature softer dual-density foam on the backrest, a re-designed pelvic well and Stretch-Air™, a breathable, water-resistant cushion fabric that is machine-washable at 60°C.

Active Reach

Tilts the seat forward by up to 45° and raises it to extend functional reach, aid in transfers and your daily activities.

Active Height

Enables you to drive up to 5 km/h while fully elevated and combines seamlessly with *Active Reach* to give you the ultimate advantage in accessing your environment.

Power recline


Smooth and safe power recline with seamless sheer reduction for any user up to 150kg.

More drivability


The new Permobil battery and Permobil VoltPro charger work together to provide you with a greater driving distance, optimized charging times and longer battery life.

Safer driving after dark

Optional LED lights front and back for improved brightness and visibility.



**Agile suspension for superior ride
comfort, traction, and stability.**




**The power recline with
seamless sheer reduction
moves between positions**

F5 Corpus[®] Specifications

Max. speed	12 km/h
Range (G24: 85 Ah batteries)	40 km*
Max. user weight	150 kg
Obstacle climbing (with run up)	65 mm (75 mm)
Min. turning diameter	1490 mm
Turning in corridor	1200 mm
Suspension	Independent on all wheels
Electronics R-net	120 A
Width incl. seat	650–790 mm
Height	970–1170 mm
Min. transport length	955 mm
Min. transport height (folded backrest)	825 mm
Weight incl. batteries 85 Ah	189 kg
Weight of batteries 85 Ah	2x23 kg
Seat height (with electric seat elevator)	450–800 mm
Seat depth	370–570 mm (by 25 mm)
Seat width	420–570 mm (by 50 mm)
Backrest width	360/410/460/510 mm
Backrest height	470, 545-670 mm (by 25 mm)
Distance between armrests	380–480/480–580 mm

Personalise your ride



*Driving conditions and choice of battery type.
Measurements are approximates and subject to change without notice.

The fully independent
Agile suspension is
especially **beneficial**
for active
outdoor driving.

Permobil

F5 Corpus®

Enjoy superior vibration reduction,
ride comfort and drive performance.
Active Reach and *Active Height*
combined seamlessly to put your
world within your reach.

The F5 Corpus is a great hybrid for
indoor or outdoor use.

Customer Service in Cyprus:



m.g.
MEDICOM SERVICES LTD

F5CORPSUDATASHEETDI211216

permobil

Tel.: 77771314, 24642070, Mob.: 99455233
e-mail: m.g.medicom@cytanet.com.cy

