



ORMESA®
• MADE IN ITALY SINCE 1980 •





BRILLO



Emotions in motion.

Birillo by Ormesa is a **walker for children** with disabilities designed to allow autonomous movement in sitting position.

Its **peculiar shape and bright colors** characterizing the **self-extinguishing upholstery**, make it a fun and stimulating aid and a valuable tool for rehabilitation activities.



PRODUCT
MARKED 

BIRILLO is a MEDICAL DEVICE THAT COMPLIES WITH THE SAFETY REQUIREMENTS defined with the REGULATION (EU) 2017/745

Modeled on you

Available in 4 sizes

Birillo is an extremely helpful aid for **the exploration of the surrounding environment**, which instantly shift from being a visual and static space to a territory of discovery. **Thanks to Birillo, the child can enrich his wealth of experiences and increase his cognitive and relational repertoire, increasing body and space awareness.**

With Birillo, safety and fun go hand in hand. The triangular shape of the seat gives stability to the pelvis and trunk, favoring a functional position to the use of the lower and upper limbs. Thanks to these features, it **stimulates the correct posture of the child without inhibiting the rotation of the trunk**, favoring the management of movement through the arms. The correct posture prevents deformities and favors the generation of motor schemes, essential for a **rich and meaningful interaction with the physical and social environment.**



Birillo 1



Birillo 2



Birillo 3



Birillo 4

Technical features and adjustments



Frame

steel and painted with non-toxic epoxy varnish.
Adjustable in height for a perfect usability and with protection balls at the upper end of the rods.



Abduction block

adjustable in depth.





Triangular seat

to contain the pelvis (seat/backrest with 90° angle). It is padded and covered with plastic material, resisting to wear and tear and cleaning.



Wheels

fully-swivelling, rolling on ball-bearings. They are protected by elastic wall-save bumpers.



Handlebar

adjustable in height and depth.

Components



860 Concave headrest
Adjustable in height.



824 Transparent tray
With circular recess. Adjustable in depth.







812 Directional locks
Help the direction of the walking frame.




Birillo sizes and weight


	size 1	size 2	size 3	size 4
A: backrest height	18 cm ~ 7 in	19 cm ~ 7,4 in	21 cm ~ 8,2 in	26 cm ~ 10,2 in
B: length encumbrance	59 cm ~ 23,2 in	65 cm ~ 25,5 in	73 cm ~ 28,7 in	78 cm ~ 30,7 in
C: overall width	50 cm ~ 19,6 in	51 cm ~ 20 in	59 cm ~ 23,2 in	60 cm ~ 23,6 in
D: height ground to seat	from 20 to 28 cm from 7,8 to 11 in	from 23 to 31 cm from 9 to 12,2 in	from 23 to 35 cm from 9 to 13,7 in	from 22 to 40 cm from 8,6 to 15,7 in
E: height ground to handlebar	from 34 to 48 cm from 13,3 to 18,8 in	from 40 to 56 cm from 15,7 to 22 in	from 40 to 65 cm from 15,7 to 25,5 in	from 43 to 74 cm from 16,9 to 29,1 in
F: seat depth	18 cm ~ 7 in	20 cm ~ 7,8 in	23 cm ~ 9 in	30 cm ~ 11,8 in
G: seat width	22 cm ~ 8,6 in	26 cm ~ 10,2 in	30 cm ~ 11,8 in	38 cm ~ 14,9 in
H: user height	from 70 to 95 cm from 27,5 to 37,4 in	from 85 to 105 cm from 33,4 to 41,3 in	from 95 to 120 cm from 37,4 to 47,2 in	from 110 to 140 cm from 43,3 to 55,1 in
Turning width	67 cm ~ 26,3 in	73 cm ~ 28,7 in	81 cm ~ 31,8 in	89 cm ~ 35 in
Total weight	8,5 kg ~ 18,7 lb	10 kg ~ 22 lb	11 kg ~ 24,2 lb	12,8 kg ~ 28,2 lb
Maximum load	30 kg ~ 66,1 lb	40 kg ~ 88,1 lb	50 kg ~ 110,2 lb	60 kg ~ 132,2 lb





The
pleasure of
discovery





Ergonomics
and
safety




m.g.

MEDICOM SERVICES LTD

Τηλ.: 77771314, 24642070, Κινητό: 99455233

www.medicomstore.com

Διεύθυνση: Ηνωμένων Πολιτειών Αμερικής 7, 6043, Λάρνακα Κύπρος



 **complies with MDR (EU) 2017/745**

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV
[ISO 13485](https://www.iso.org/standard/54554.html)

Ormesa srl - Via delle Industrie 6/8 | 06034 Foligno (PG) - ITALY
Tel. +39 0742 22927 | FAX +39 0742 22637 | info@ormesa.com

www.ormesa.com

