

Permobil

Panthera S3 Swing



permobil

Adapts to the active life you lead

The Panthera S3 Swing offers the same comfort and low transport weight as the Panthera S3 series, the S3 Swing provides the added flexibility of a swing-away foot rest for easier transfers. Also available in S3 Swing Short and Long option to suit smaller or taller builds.

Key benefits:

- Made to maneuver
- Low weight and robust
- Adjusts for comfort
- Integrated anti-tip device
- Easy to maintain



Specifications

	S3 Swing	S3 Swing Short
Seat width	36 - 39 - 42 - 45 - 50 cm	33 - 36 - 39 cm
Seat depth	40 cm	35 cm
Seat height	Front 45 cm (24" rear wheels) Back 43 cm (24" rear wheels)	Front 45 cm (24" rear wheels) Back 43 cm (24" rear wheels)
Seat angle	4°	4°
Total length	78 - 90 cm	88 cm
Total width	Seat width +21 cm	Seat width +21 cm
Front frame angle	105°	105°
Backrest height	20 - 45 cm	20 - 45 cm
Backrest angle	13° - (-8°)	13° - (-8°)
Footplate distance	38 - 46 cm	38 - 46 cm
Camber	2.2°, 4° (seat width 36 - 45 cm) 8° (seat width 36 - 39 cm)	2.2°, 4°, 6 (seat width 33 cm), 8°
Weight without rear wheels	4098 - 4898 grams	3958 - 4078 grams
Wheelchair weight	9029 - 9957 grams	8850 - 9006 grams
Maximum load	100 kg, 150 kg (seat width 45 - 50 cm)	100 kg
Material (chassis)	Chrome molybdenum tubing	Chrome molybdenum tubing
Material (upholstery)	Polyurethane covered polyester	Polyurethane covered polyester

Personalise your ride*

ANODIC BLACK

PEARL WHITE

GREY METALLIC

BLUE CLASSIC

PINE GREEN

Permobil reserves the right to change specifications and design without prior notice. *More color options available.

Actual product colors may vary.