

Permobil

Panthera U3 Light



permobil

Optimal ergonomics to enable your lifestyle

A low-weight, maneuverable manual wheelchair, the Panthera U3 Light offers stability and comfortable posture support. Designed for optimal ergonomics when seated or driving, this wheelchair can also be easily lifted into a car thanks to its balanced, effortless-grip chassis. Available as U3 Light and U3 Light L.

Key benefits:

- Ultra-light weight
- Extremely maneuverable
- Easy to transport
- Adaptable to individual needs
- Easy to maintain



Specifications

| | U3 Light | U3 Light Y-front |
|----------------------------|---|---|
| Seat width | 33 - 36 - 39 - 42 - 45 cm | 33 - 36 - 39 - 42 - 45 cm |
| Seat depth | 35 - 46 cm | 35 - 46 cm |
| Seat height | Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels) | Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels) |
| Seat angle | 7° | 7° |
| Total length | 82 cm | 82 cm |
| Total width | Seat width +21 cm | Seat width +21 cm |
| Front frame angle | 100° | 100° |
| Backrest height | 25 - 35 cm | 25 - 35 cm |
| Backrest angle | 11.5° - (-4.5°) | 11.5° - (-4.5°) |
| Footplate distance | 35 - 43 cm | 35 - 43 cm |
| Camber | 2.2° | 2.2° |
| Weight without rear wheels | 3.4 - 4.0 kg | 3.7 - 3.9 kg |
| Wheelchair weight | 5.7 - 6.3 kg | 6.0 - 6.2 kg |
| Maximum load | 100 kg | 100 kg |
| Material (chassis) | Chrome molybdenum tubing | Chrome molybdenum tubing |
| Material (upholstery) | Polyurethane covered polyester | Polyurethane covered polyester |



Personalise your ride*

ANODIC
BLACK

PEARL
WHITE

GREY
METALLIC

BLUE
CLASSIC

PINE
GREEN

Actual product colors may vary.

U3 Light Y-front 90°

33 - 36 - 39 - 42 - 45 cm

35 - 46 cm

Front 47 cm (24" rear wheels)
Back 43 cm (24" rear wheels)

7°

73 cm

Seat width +21 cm

90°

25 - 35 cm

11.5° - (-4.5°)

35 - 43 cm

2.2°

4.0 - 4.2 kg

6.3 - 6.5 kg

100 kg

Chrome molybdenum tubing

Polyurethane covered polyester

U3 Light Y-front High 4°

33 - 36 - 39 - 42 - 45 cm

35 - 46 cm

Front 51 cm (24" rear wheels)
Back 45 cm (24" rear wheels)

4°

82 cm

Seat width +21 cm

100°

25 - 35 cm

11.5° - (-4.5°)

36 - 44 cm

2.2°

3.7 - 3.9 kg

6.0 - 6.2 kg

100 kg

Chrome molybdenum tubing

Polyurethane covered polyester

U3 Light Y-front Low

33 - 36 - 39 - 42 - 45 cm

35 - 46cm

Front 44.5 cm (w. 24" rear wheels)
Back 40.5 cm (w. 24" rear wheels)

7°

82 cm

Seat width +21 cm

100°

25 - 35 cm

11.5° - (-4.5°)

36 - 44 cm

2.2°

3.7 - 3.9 kg

6.0 - 6.2 kg

100 kg

Chrome molybdenum tubing

Polyurethane covered polyester



Enjoy
the ride

Panthera
U3 Light

For experienced manual wheelchair users who need a low-weight, maneuverable mobility device, the Panthera U3 Light offers stability and comfortable posture support.

Pair it with Mototronik* power assist device for added convenience.

*Mototronik - see picture

permobil



PANTHERAU3LIGHTDATASHEETD1221228