

Permobil

# Panthera U3 Series



permobil

# Considered design that lets you live your life your way

The Panthera U3 is a low-weight, adjustable manual wheelchair that offers a feeling of security in an ergonomically optimal driving position. Made for active users who want a slim-looking wheelchair and easy maneuverability, combined with adjustable accessories for added comfort and security.

## Key benefits:

- Light and strong
- Adjusts to individual needs
- Integrated anti-tip device
- Easy to manage



## Specifications

	U3	U3 Y-front
Seat width	33 - 36 - 39 - 42 - 45 cm	33 - 36 - 39 - 42 - 45 - 50 cm
Seat depth	35 - 46 cm	35 - 46 cm
Seat height	Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels)	Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels)
Seat angle	7°	7°
Total length	84 cm	84 cm
Total width	Seat width +21 cm	Seat width +21 cm
Front frame angle	100°	100°
Backrest height	20 - 45 cm	20 - 45 cm
Backrest angle	13° - (-8°)	13° - (-8°)
Footplate distance	35 - 43 cm	35 - 43 cm
Camber	2.2°, 4°, 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)	2.2°, 4°, 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)
Weight without rear wheels	4.5 - 5.0 kg	4.5 - 5.4 kg
Wheelchair weight	8.0 - 8.5 kg	8.0 - 9.3 kg
Maximum load	100 kg, 150 kg (seat width 45 cm)	100 kg, 150 kg (seat width 45 - 50 cm)
Material (chassis)	Chrome molybdenum tubing	Chrome molybdenum tubing
Material (upholstery)	Polyurethane covered polyester	Polyurethane covered polyester



Personalise your ride\*

- ANODIC BLACK
- PEARL WHITE
- GREY METALLIC
- BLUE CLASSIC
- PINE GREEN

Actual product colors may vary.

U3 Y-front Short	U3 Y-front Low
33 - 36 - 39 - 42 - 45 - 50 cm	33 - 36 - 39 - 42 - 45 - 50 cm
35 - 46 cm	35 - 46 cm
Front 47 cm (24" rear wheels) Back 43 cm (24" rear wheels)	Front 44.5 cm (24" rear wheels) Back 40.5 cm (24" rear wheels)
7°	7°
79 cm	84 cm
Seat width +21 cm	Seat width +21 cm
100°	100°
20 - 45 cm	20 - 45 cm
13° - (-8°)	13° - (-8°)
35 - 43 cm	35 - 43 cm
2.2°, 4°, 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)	2.2°, 4°, 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)
4.5 - 5.4 kg	4.4 - 5.3 kg
8.0 - 9.3 kg	8.0 - 9.3 kg
100 kg, 150 kg (seat width 45 - 50 cm)	100 kg, 150 kg (seat width 45 - 50 cm)
Chrome molybdenum tubing	Chrome molybdenum tubing
Polyurethane covered polyester	Polyurethane covered polyester

# Innovating for individuals