

Permobil

# Progeo Noir 2.0

Performance and design in perfect harmony



Carbon



Rigid



Minimized  
weight



Made-to-  
measure

permobil

# The perfect fit

Progeo Noir 2.0 combines the benefits of carbon with the superior manoeuvrability of a made to measure wheelchair. This together gives the user an optimal balanced ride allowing for energy saving throughout the day.



## Progeo Noir 2.0

Seat width	33 - 36 - 39 - 42 - 45 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm
Seat height	Made-to-measure
Footplate distance	Made-to-measure (standard) Adjustable (on request)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	85° - 90° - 95°
Total length	85 cm approx. (standard configuration)
Total width	Seat width + 17 cm (with 0° camber)
Backrest angle	From 74° to 94° to the seat (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 4,0 kg approx.
Wheelchair weight	Starting from 6,6 kg approx.
Maximum load	125 kg

CARBON  
ESSENCE

LIME  
ESSENCE

RED  
ESSENCE

BLUE  
ESSENCE

OPAQUE  
ESSENCE

Personalise your ride

Permobil reserves the right to change specifications and design without prior notice.

Actual product colors may vary.

permobil

Customer Service in Cyprus:  
 m.g.  
MEDICOM SERVICES LTD

Tel.: 77771314, 24642070, Mob.: 99455233  
e-mail: m.g.medicom@cytanet.com.cy



PROGEOINOIR2.0DATASHEETD1210909